



Patient Name:		Phone:	
Date of birth:		Diagnosis/ICD-10 Code(s):	
Selecting a profile results in the allergens below it being tested. Allergens listed w/reflex that are positive (>0.10) will result in the <i>components</i> below it being tested. Individual allergens can be selected instead of or in addition to profiles. Reflex can be overridden by request in the comment section.			
Respiratory Allergy Profile Reg 6	Food Allergy Profile	Child/Pediatric Allergy Profile	Miscellaneous Allergens
Order Code: ALLERGENRESP	Order Code: ALLERGENFOOD	Order Code: ALLERGENCHILD	Gluten
☐ D. Pteronyssinus	☐ Peanut w/reflex	☐ D. Pteronyssinus	☐ Sweet Vernal Grass
☐ D. Farinae	☐ <i>Ara h 1</i>	☐ D. Farinae	☐ Meadow Fescue
☐ Cat Dander w/reflex	☐ <i>Ara h 2</i>	☐ Cat Dander w/reflex*	☐ Rye Grass
☐ <i>Fel d 1</i>	☐ <i>Ara h 3</i>	☐ Dog Dander w/reflex*	☐ Meadow, Kentucky Blue Grass
☐ <i>Fel d 2</i>	☐ <i>Ara h 6</i>	☐ Mouse	☐ Johnson Grass
☐ <i>Fel d 4</i>	☐ <i>Ara h 8</i>	☐ Cockroach, German	☐ Bahia Grass
☐ <i>Fel d 7</i>	☐ <i>Ara h 9</i>	☐ Penicillium Chrysogenum	☐ Latex
☐ Dog Dander w/reflex	☐ MUXF3 CCD – Bromelain	☐ Cladosporium Herbarum	☐ Candida Albicans
☐ <i>Can f 1</i>	☐ <i>Bet v 2 (profilin)</i>	☐ Cow 's Milk w/reflex*	☐ Helminthosporia
☐ <i>Can f 2</i>	☐ Hazelnut w/reflex	☐ Egg White w/reflex*	☐ Fusarium Proli
☐ <i>Can f 3</i>	☐ <i>Cor a 1</i>	☐ Peanut w/reflex*	☐ Stemphylium Herbarum
☐ <i>Can f 4</i>	☐ <i>Cor a 8</i>	☐ Wheat w/reflex*	☐ Rhizopus Nigricans
☐ <i>Can f 5</i>	☐ <i>Cor a 9</i>	☐ <i>Gliadin w/reflex*</i>	☐ Epicoccum Purpurascens
☐ <i>Can f 6</i>	☐ <i>Cor a 14</i>	☐ Soybean	☐ Curvularia Lunata
☐ Mouse	☐ MUXF3 CCD – Bromelain	☐ Cod Fish	☐ Acremonium Kiliense
☐ Cockroach, German	☐ <i>Bet v 2 (profilin)</i>	☐ Shrimp	☐ Trichophyton Rubrum
☐ Penicillium Chrysogenum	☐ Sesame Seed w/reflex	☐ Total IgE	☐ Phoma
☐ Cladosporium Herbarum	☐ <i>Ses i 1</i>		☐ Plantain
☐ Aspergillus Fumigatus	☐ Brazil Nut w/reflex	Infant Food Profile	☐ Goosefoot, Lambs Quarters
☐ Alternaria Alternata	☐ <i>Ber e 1</i>	Order Code: INFANTFOOD	☐ Sheep Sorrel
☐ Timothy Grass	☐ Cashew w/reflex	☐ Cow 's Milk w/reflex*	☐ Red Cedar
☐ Bermuda Grass	☐ <i>Ana o 3</i>	☐ Egg White w/reflex*	☐ Red Mulberry
☐ Box-elder/Maple	☐ Walnut w/reflex	☐ Peanut w/reflex*	☐ Russian Thistle
☐ Walnut	☐ <i>Jug r 1</i>	☐ Wheat w/reflex*	☐ Cocklebur
☐ Pecan/Hickory	☐ <i>Jug r 3</i>	☐ <i>Gliadin w/reflex*</i>	☐ Sycamore/London Plane
☐ Common Silver Birch	☐ Cow 's Milk w/reflex	☐ Soybean	☐ Sweet Gum
☐ Mountain Juniper	☐ <i>α-lactalbumin</i>		☐ Cottonwood
☐ Oak	☐ <i>β-lactoglobulin</i>	*Reflex components listed under	☐ Ash/White
☐ Mulberry	☐ <i>Casein</i>	the Respiratory and Food Profiles	☐ Pine
☐ Elm	☐ Egg White w/reflex	Others/Comments:	
☐ Common Ragweed	☐ <i>Ovomucoid</i>		
☐ Common Pigweed	☐ <i>Ovalbumin</i>		
☐ Rough Marshelder	☐ Wheat w/reflex		
☐ Total IgE	☐ <i>Tri a 14</i>		
	☐ MUXF3 CCD (Bromelain)		
Alpha Gal Profile	☐ <i>Phl p 12 (profilin)</i>		
Order Code: ALPHAGAL	☐ <i>Gliadin w/reflex</i>		
☐ alpha-Gal	☐ <i>Tri a 19 (omega-5-gliadin)</i>		
☐ Beef	☐ Soybean		
☐ Lamb	☐ Almond		
☐ Pork	☐ Cod Fish		
☐ Total IgE	☐ Salmon		
	☐ Tuna		
	☐ Shrimp		
	☐ Scallop		
	☐ Total IgE		
Whole allergen-CPT 86003		Specimen Type: Serum is the sample of choice. Collect 2 gold tops. Sample should be stored at 2-8 degrees C if assayed within one week. Otherwise specimens should be frozen at -20 degrees C. No specific patient preparation is required.	
Component -CPT 86008			
Total IgE-CPT 82785			
Provider's Signature:		Date:	